

Kaye Products, Inc.

ADAPTIVE EQUIPMENT AND REHABILITATION PRODUCTS
535 DIMMOCKS MILL ROAD HILLSBOROUGH, NC 27278 USA
TEL. (919) 732-6444 FAX (800) 685-5293

KAYE ADJUSTABLE BENCHES

(effective 4/11/2018)

These are the user instructions for the Kaye Adjustable Benches. All people who set-up, adjust, and use the bench should first read these entire instructions.

Please store these instructions for future reference.

USE

These benches are designed as a general purpose sitting device. They can serve as a positioning tool for users who do not need upper body support as they sit, but who may need a stable, yet adjustable seat to help position their lower body while they sit.

BENCH DESIGNS

The S1A, S2A, S3A, S4A, and S5A have wider cushions. The extra width provides added stability and room for two users.

The TSS and TSL, with their narrower cushions, are designed to fit under small desks and in tight spaces.

The S3A and S4A fold flat for easy transportation and storage.

BENCH SIZES

The model number is located on the warning label. For older benches, you can determine the model by comparing the bench's dimensions to the following chart. When referring to the chart:

A = The cushion size (depth times width).

Inches are on the top row, centimeters on the bottom row.

B = The height adjustment from the floor to the top of the seat. Inches are on the top row, centimeters on the bottom row.

C = Whether the bench folds.

Model	A (ins./cms)	B (ins./cms)	C (folds)
S1A	10 1/2 x 26, 27 x 66	7 - 11, 18 - 28	no
S2A	14 x 31, 36 x 79	10 1/2 - 17, 27 - 43	no
S3A	10 1/2 x 26, 27 x 66	7 - 11, 18 - 28	yes
S4A	14 x 31, 39 x 79	10 1/2 - 17, 27 - 43	yes
S5A	14 x 31, 39 x 79	16 - 23, 41 - 58	no
TSS	10 1/2 x 12, 27 x 30	7 - 11, 18 - 28	no
TSL	14 x 17 1/2, 39 x 44	10 1/2 - 17, 27 - 43	no

FOLDING AND UNFOLDING THE S3A AND S4A

To unfold the S3A or S4A, rotate the two legs down until the latch push-pins protrude from the corresponding slot in the telescoping tubes.

To fold the bench, depress the latch push-pins and rotate the leg assemblies in until they rest against the underside of the bench.

ASSEMBLY


Some benches are shipped with the legs removed and wrapped separately. Install the legs as follows:

1. Remove the black knobs and washers from each end of the bench.
2. Install the legs with the labels facing outward and the rubber feet pointing downward
3. Place the washers and then the knobs onto the bolts and tighten the knobs to secure the leg in place.

TILT AND HEIGHT ADJUSTMENT

To tilt and adjust the height of the bench, adjust each leg as follows:

1. Loosen the black knobs.
2. Slide the leg sideways until the bolts are in the vertical slots.
3. Move the leg up or down.
4. To tilt the leg, move the front edge of the leg up or down relative to the rear as shown:



Rear

Front

5. Slide the leg sideways to seat the bolts fully in the ends of the horizontal slots.
6. Tighten the knobs.

ACCESSORIES

A posture system is available for each of the non-folding benches. That system includes a pelvic support, a femoral stabilizer, and an abductor. When added to the bench, this system provides added stability to the pelvis and spine.

The accessories include their own user instructions.

Please call Kaye Products or consult their catalog for more details.

MAINTENANCE AND CARE

Inspect the unit regularly. Tighten the nuts and bolts as necessary.

If the unit needs service or spare parts, contact Kaye Products, Inc. or the distributor from whom you purchased the item. If a problem is discovered that may impact the unit's function, immediately cease use and contact Kaye Products, Inc.

Do not expose the unit to rain or submerge it in water.

Use a slightly damp cloth to clean the unit, but wipe the dampness from the wood after cleaning it.

Avoid any undue stress to the bench while using, storing, or transporting it.

LIMITED WARRANTY

If an item proves defective within two years of the original purchase, we will provide replacement parts in order to correct that defect. Normal wear and tear is not covered by the warranty.

Kaye Products, Inc. makes no other warranty, expressed or implied, and does not warrant the product as being fit for a particular purpose. The purchaser, owner, and user assume all risk of personal and property injury due to the use of the equipment.

CAUTIONS

1. Do not fold or adjust the bench while it is in use.
2. Each time you unfold the S3A or S4A, ensure that the two push-pins fully extend from their corresponding slots.
3. Do not use if there are broken or missing parts.
4. The maximum weight limit for all models is 230 pounds or 104 kilograms.
5. Do not alter the bench or use it in any way other than described in these instructions.
6. Before use, always ensure that the black knobs are tightened and the legs are adjusted so that the bolts rest fully in the horizontal slots.
7. Do not leave a child unattended.
8. Do not tilt the leg so that one knob is positioned more than two slots higher or lower than the other.
9. Do not stand on the bench.

QUESTIONS?

If you have any questions, comments, or requests, please contact Kaye Products, Inc.



Our authorized representative for the European Medical Device Directives is Advena, Ltd., PO Box 30, Leominster, Herefordshire, HR6 0ZQ, UK.