

Some styles of Kaye Posture Control Walkers can be converted to an anterior configuration.

Walker Style:	Conversion Kit
W1/2B	W1/2BC-ANT
W1/2BR	W1/2BRC-ANT
W1B	W1BC-ANT
W1BR	W1BRC-ANT
W2B	No kit needed on these styles. Simply switch front and rear legs.
W2BR	
W3B	
W3BR	
W4B	
W4BR	
W5C	

It is NOT possible to make a swivel wheeled 'B' and 'C' Frame walker work as an anterior walker.

All 'H' Frame Walkers can be made anterior simply by switching the legs. This includes swivel-wheeled 'H' walkers.

Model No.	Y1FS	Y2FS	Y3FS	Y4FS
Weight limits	85 lb. (39 kg.)	85 lb. (39 kg.)	130 lb. (60 kg.)	180 lb. (82 kg.)
Height to top of forearm cushion (floor to elbow)	19" to 26"	21½" to 30"	28" to 39½"	37" to 45"
Chest cushions	5" w x 6" h	6½" w x 8" h	7½" w x 9¼" h	7½" w x 9¼" h
Greatest outside dimensions	25" l x 24" w	32" l x 26" w	37" l x 26" w	45" l x 26" w
Distance between forearm cushions	6" to 10"	6½" to 12½"	10" to 17"	12½" to 19"
Distance between trunk support cushions	6" to 10"	8" to 12"	9" to 15"	10" to 15½"
Distance between forearm cushion and handgrip	1½" to 5"	1½" to 5½"	1½" to 5½"	3" to 7"
Accessories				
Swivel Limiters	WSL-1	WSL-1	WSL-1	WSL-1
Variable Resistance Rear Wheels	VRY1	VRY2	VRY3	VRY4
Soft Sling Support	SSS1	SSS1	SSS2	SSS2

Model No.	Y1S	Y2S	Y3S
Weight limits	85 lb. (39 kg.)	85 lb. (39 kg.)	130 lb. (60 kg.)
Height to top of forearm cushion (floor to elbow)	19" to 26"	22" to 36"	30" to 47"
Chest cushions	5" w x 6" h	6½" w x 8" h	7½" w x 9½" h
Greatest outside dimensions	25" l x 24" w	32" l x 26" w	37" l x 26" w
Distance between lateral support cushions	6" to 11"	9" to 14"	10" to 15"
Accessories			
Swivel Limiters	WSL-1	WSL-1	WSL-1
Variable Resistance Rear Wheels	VRY1	VRY2	VRY3
Soft Sling Support	SSS1	SSS1	SSS2

Walker Accessories

- Swivel Limiters (Page 6)
- Variable Resistance Rear Wheels (Page 7)
- Soft Sling Support (Page 18)



Y4FS