

How to Choose your Kaye Posture Control Walker

1. Choose the type of frame that you want:

- Posture Control 'B' Frame Walker (pg 4)
- Posture Control 'R' Wide Frame Walker (pg 4)
- Posture Control 'C' Frame Walker (pg 6)
- Posture Rest 'H' Frame Walker (pg 6)
- My Red Walker (pg 5)

2. Choose the size walker that you need, sizes are indicated by numerals: (Chart No. 1 – below)

3. Choose your Walker Configuration:

- Two Wheeled Walkers: with and without seat (Chart No. 2 – page 21)
- Four Wheeled Walkers:
 - Non-Swivel Front Wheels with and without seat (Chart No. 4 – page 22)
 - Swivel Front Wheels with and without seat (Chart No. 5 – page 23)

4. Choose your accessories: (Charts No. 2 – No. 5)

Chart 1: Posterior Walker Frame Measurements

Model No.	Height to top of handle (floor to mid-buttocks)	Distance inside handles (depth/width)	Weight Limits (lbs/kg)
Posture Control 'B' and 'C' Frame Walkers			
W1/2B*	14.5" to 20"	9" / 13.5"	60 / 27
W1B	16.5" to 22"	9" / 13.5"	60 / 27
W2B	19" to 25"	9.5" / 13.5"	85 / 39
W3B	23" to 30.5"	10.5" / 15"	130 / 60
W4B	28.5" to 36"	13.5" / 16"	210 / 95
W5C	35.5" to 41"	15" / 18.5"	250 / 113
*Height to Top of Handle on W1/2B with rear wheels is 14½" - 18½"			
Posture Control 'R' Frame Wide Walkers			
R2B	19" to 25"	10.5" / 15.5"	85 / 39
R3B	23" to 30."	12.5" / 17.5"	130 / 60
R4B	28.5" to 36"	13.5" / 19"	210 / 95
Posture Rest 'H' Frame Walkers			
W1/2BH	14.5" to 20"	9" / 13.5"	60 / 27
W1BH	16.5" to 22"	9" / 13.5"	60 / 27
W2H	20" to 26.5"	11" / 16.25"	85 / 39
W3H	25.5" to 32"	12" / 19"	130 / 60
W4H	31" to 37"	14" / 21"	180 / 82
W5H	36" to 41"	13.5" / 20"	250 / 113
W5CH	35.5" to 41"	15" / 18.5"	250 / 113
My Red Walker / Walker Wagon			
RFA	13.5" to 18.5"	9.5" / 13.5"	60 / 27
RFB	17" to 25"	9.5" / 13.5"	85 / 39
WW1	15" to 25"	14.75"w	
WW2	21" to 31"	14.75"w	

Choose the correct walker and accessories quickly and easily with our walker selector app, available online at kayeproducts.com/web-app/.



Scan code to access the online app.



Chart 2: Two Wheeled Walkers

(Non-swivel front, crutch tips on rear)

Model No.	Add-A-Seat (p. 7)	Leg Abductor (p. 14)	Suspension Conversion Kit (p. 14)	All Terrain Wheels (p. 8)	Silent Rear Wheels (p. 9)	Variable Resistance Wheels (p. 9)	5" Ratchet Rear Wheels (p. 9)	Forearm Supports (p. 13)	Extensor Assist Pad (p. 12)	Pelvic Stabilizer Pads (p. 12)
Non-Swivel Front Wheels (R), Rear Crutch Tips										
W1/2B	W1S	LA1	SC1	N/A	RR1/2B	VR1/2B	N/A	FSS	EA2	PS2H
W1B	W1S	LA1	SC1	W1BAC	RR1B	VR1B	W1BAR	FSS	EA2	PS2H
W2B	W2S	LA2	SC2	W2BAC	RR2B	VR2B	W2BAR	FSS	EA2	PS2H
R2B	N/A	N/A	SC2R	W2BAC	RR2B	VR2B	W2BAR	FSS	EA3	PS3H
W3B	W3S	LA3	SC3	W3BAC	RR3B	VR3B	W3BAR	FSS	EA3	PS3H
R3B	N/A	N/A	SC3R	W3BAC	RR3B	VR3B	W3BAR	FSS	EA3	PS4H
W4B	W4S	LA4	SC4	W4BAC	RR4B	VR4B	W4BAR	FSM	EA3	PS3H
R4B	N/A	N/A	SC4R	W4BAC	RR4B	VR4B	W4BAR	FSM	EA3	PS4H
W5C	5CH	N/A	N/A	W4BAC	RR4B	VR5C	W4BAR	FSM	EA3	PS4H
Non-Swivel Front Wheels (R), Rear Crutch Tips, Seat (H)										
W1/2BH	included	LA1	SC1	N/A	RR1/2B	VR1/2B	N/A	FSS	EA2	PS2H
W1BH	included	LA1	SC1	W1BAC	RR1B	VR1B	W1BAR	FSS	EA2	PS2H
W2H	included	N/A	N/A	W2HAC	RR2H	VR2H	W2HAR	FS2H	EA2	PS2H
W3H	included	N/A	N/A	W3HAC	RR3H	VR3H	W3HAR	FS3H	EA3	PS3H
W4H	included	N/A	N/A	W4HAC	RR4H	VR4H	W4HAR	FS4H	EA3	PS4H
W5H	included	N/A	N/A	W5HAC	RR5H	VR5H	W5HAR	N/A	EA3	PS4H
W5CH	included	N/A	N/A	W4BAC	RR4B	VR5C	W4BAR	FSM	EA3	PS4H

Chart 3: My Red Walker

(only available with swivel front, ratchet rear)

Model No.	Add-A-Seat (p. 7)	Leg Abductor (p. 14)	Suspension Conversion System (p. 14)	All Terrain Wheels (p. 8)	Silent Rear Wheels (p. 9)	Variable Resistance Wheels (p. 9)	5" Ratchet Rear Wheels (p. 9)	Forearm Supports (p. 13)	Extensor Assist Pad (p. 12)	Pelvic Stabilizer Pads (p. 12)
RFA	N/A	N/A	N/A	N/A	N/A	N/A	N/A	FSA	EA2	PS2H
RFB	N/A	N/A	N/A	N/A	N/A	N/A	N/A	FSA	EA2	PS2H

Chart 4: Four Wheeled Walkers Non-Swivel Front Wheels

Model No.	Add-A-Seat (p. 7)	Leg Abductor (p. 14)	Suspension Conversion Kit (p. 14)	All Terrain Wheels (p. 8)	Silent Rear Wheels (p. 9)	Variable Resistance Wheels (p. 9)	5" Ratchet Rear Wheels (p. 9)	Forearm Supports (p. 13)	Extensor Assist Pad (p. 12)	Pelvic Stabilizer Pads (p. 12)
Non-Swivel Front Wheels (R), One-Way Ratchet Rear Wheels										
W1/2BR	W1/2S	LA1	SC1	N/A	RR1/2B	VR1/2B	N/A	FSS	EA2	PS2H
W1BR	W1S	LA1	SC1	W1BAC	RR1B	VR1B	W1BAR	FSS	EA2	PS2H
W2BR	W2S	LA2	SC2	W2BAC	RR2B	VR2B	W2BAR	FSS	EA2	PS2H
R2BR	N/A	N/A	SC2R	W2BAC	RR2B	VR2B	W2BAR	FSS	EA3	PS3H
W3BR	W3S	LA3	SC3	W3BAC	RR3B	VR3B	W3BAR	FSS	EA3	PS3H
R3BR	N/A	N/A	SC3R	W3BAC	RR3B	VR3B	W3BAR	FSS	EA3	PS4H
W4BR	W4S	LA4	SC4	W4BAC	RR4B	VR4B	W4BAR	FSM	EA3	PS3H
R4BR	N/A	N/A	SC4R	W4BAC	RR4B	VR4B	W4BAR	FSM	EA3	PS4H
W5CR	5CH	N/A	N/A	W4BAC	RR4B	VR5C	W4BAR	FSM	EA3	PS4H
Non-Swivel Front Wheels (R), One-Way Ratchet Rear Wheels, Seat (H)										
W1/2BHR	included	LA1	SC1	N/A	RR1/2B	VR1/2B	N/A	FSS	EA2	PS2H
W1BHR	included	LA1	SC1	W1BAC	RR1B	VR1B	W1BAR	FSS	EA2	PS2H
W2HR	included	N/A	N/A	W2HAC	RR2H	VR2H	W2HAR	FS2H	EA2	PS2H
W3HR	included	N/A	N/A	W3HAC	RR3H	VR3H	W3HAR	FS3H	EA3	PS3H
W4HR	included	N/A	N/A	W4HAC	RR4H	VR4H	W4HAR	FS4H	EA3	PS4H
W5HR	included	N/A	N/A	W5HAC	RR5H	VR5H	W5HAR	N/A	EA3	PS4H
W5CHR	included	N/A	N/A	W4BAC	RR4B	VR5C	W4BAR	FSM	EA3	PS4H
Non-Swivel Front Wheels (R), Silent Rear Wheels (X)										
W1/2BRX	W1/2S	LA1	SC1	N/A	included	VR1/2B	N/A	FSS	EA2	PS2H
W1BRX	W1S	LA1	SC1	W1BAC	included	VR1B	W1BAR	FSS	EA2	PS2H
W2BRX	W2S	LA2	SC2	W2BAC	included	VR2B	W2BAR	FSS	EA2	PS2H
R2BRX	N/A	N/A	SC2R	W2BAC	included	VR2B	W2BAR	FSS	EA3	PS3H
W3BRX	W3S	LA3	SC3	W3BAC	included	VR3B	W3BAR	FSS	EA3	PS3H
R3BRX	N/A	N/A	SC3R	W3BAC	included	VR3B	W3BAR	FSS	EA3	PS4H
W4BRX	W4S	LA4	SC4	W4BAC	included	VR4B	W4BAR	FSM	EA3	PS3H
R4BRX	N/A	N/A	SC4R	W4BAC	included	VR4B	W4BAR	FSM	EA3	PS4H
W5CRX	5CH	N/A	N/A	W4BAC	included	VR5C	W4BAR	FSM	EA3	PS4H
Non-Swivel Front Wheels (R), Silent Rear Wheels (X), Seat (H)										
W1/2BHRX	included	LA1	SC1	N/A	included	VR1/2B	N/A	FSS	EA2	PS2H
W1BHRX	included	LA1	SC1	W1BAC	included	VR1B	W1BAR	FSS	EA2	PS2H
W2HRX	included	N/A	N/A	W2HAC	included	VR2H	W2HAR	FS2H	EA2	PS2H
W3HRX	included	N/A	N/A	W3HAC	included	VR3H	W3HAR	FS3H	EA3	PS3H
W4HRX	included	N/A	N/A	W4HAC	included	VR4H	W4HAR	FS4H	EA3	PS4H
W5HRX	included	N/A	N/A	W5HAC	included	VR5H	W5HAR	N/A	EA3	PS4H
W5CHRX	included	N/A	N/A	W4BAC	included	VR5C	W4BAR	FSM	EA3	PS4H

Chart 5: Four Wheeled Walkers Swivel Front Wheels

Model No.	Add-A-Seat (p. 7)	Leg Abductor (p. 14)	Suspension Conversion Kit (p. 14)	All Terrain Wheels (p. 8)	Silent Rear Wheels (p. 9)	Variable Resistance Wheels (p. 9)	5" Ratchet Rear Wheels (p. 9)	Forearm Supports (p. 13)	Extensor Assist Pad (p. 12)	Pelvic Stabilizer Pads (p. 12)
Swivel Front Wheels (S), One-Way Ratchet Rear Wheels										
W1/2BS	W1/2S	LA1	SC1	N/A	RR1/2B	VR1/2B	N/A	FSS	EA2	PS2H
W1BS	W1S	LA1	SC1	W1BAC	RR1B	VR1B	W1BAR	FSS	EA2	PS2H
W2BS	W2S	LA2	SC2	W2BAC	RR2B	VR2B	W2BAR	FSS	EA2	PS2H
R2BS	N/A	N/A	SC2R	W2BAC	RR2B	VR2B	W2BAR	FSS	EA3	PS3H
W3BS	W3S	LA3	SC3	W3BAC	RR3B	VR3B	W3BAR	FSS	EA3	PS3H
R3BS	N/A	N/A	SC3R	W3BAC	RR3B	VR3B	W3BAR	FSS	EA3	PS4H
W4BS	W4S	LA4	SC4	W4BAC	RR4B	VR4B	W4BAR	FSM	EA3	PS3H
R4BS	N/A	N/A	SC4R	W4BAC	RR4B	VR4B	W4BAR	FSM	EA3	PS4H
W5CS	5CH	N/A	N/A	W4BAC	RR4B	VR5C	W4BAR	FSM	EA3	PS4H
Swivel Front Wheels (S), One-Way Ratchet Rear Wheels, Seat (H)										
W1/2BHHS	included	LA1	SC1	N/A	RR1/2B	VR1/2B	N/A	FSS	EA2	PS2H
W1BHHS	included	LA1	SC1	W1BAC	RR1B	VR1B	W1BAR	FSS	EA2	PS2H
W2HHS	included	N/A	N/A	W2HAC	RR2H	VR2H	W2HAR	FS2H	EA2	PS2H
W3HHS	included	N/A	N/A	W3HAC	RR3H	VR3H	W3HAR	FS3H	EA3	PS3H
W4HHS	included	N/A	N/A	W4HAC	RR4H	VR4H	W4HAR	FS4H	EA3	PS4H
W5HHS	included	N/A	N/A	W5HAC	RR5H	VR5H	W5HAR	N/A	EA3	PS4H
W5CHHS	included	N/A	N/A	W4BAC	RR4B	VR5C	W4BAR	FSM	EA3	PS4H
Swivel Front Wheels (S), Silent Rear Wheels (X)										
W1/2BSX	W1/2S	LA1	SC1	N/A	included	VR1/2B	N/A	FSS	EA2	PS2H
W1BSX	W1S	LA1	SC1	W1BAC	included	VR1B	W1BAR	FSS	EA2	PS2H
W2BSX	W2S	LA2	SC2	W2BAC	included	VR2B	W2BAR	FSS	EA2	PS2H
R2BSX	N/A	N/A	SC2R	W2BAC	included	VR2B	W2BAR	FSS	EA3	PS3H
W3BSX	W3S	LA3	SC3	W3BAC	included	VR3B	W3BAR	FSS	EA3	PS3H
R3BSX	N/A	N/A	SC3R	W3BAC	included	VR3B	W3BAR	FSS	EA3	PS4H
W4BSX	W4S	LA4	SC4	W4BAC	included	VR4B	W4BAR	FSM	EA3	PS3H
R4BSX	N/A	N/A	SC4R	W4BAC	included	VR4B	W4BAR	FSM	EA3	PS4H
W5CSX	5CH	N/A	N/A	W4BAC	included	VR5C	W4BAR	FSM	EA3	PS4H
Swivel Front Wheels (S), Silent Rear Wheels (X), Seat (H)										
W1/2BHHSX	included	LA1	SC1	N/A	included	VR1/2B	N/A	FSS	EA2	PS2H
W1BHHSX	included	LA1	SC1	W1BAC	included	VR1B	W1BAR	FSS	EA2	PS2H
W2HHSX	included	N/A	N/A	W2HAC	included	VR2H	W2HAR	FS2H	EA2	PS2H
W3HHSX	included	N/A	N/A	W3HAC	included	VR3H	W3HAR	FS3H	EA3	PS3H
W4HHSX	included	N/A	N/A	W4HAC	included	VR4H	W4HAR	FS4H	EA3	PS4H
W5HHSX	included	N/A	N/A	W5HAC	included	VR5H	W5HAR	N/A	EA3	PS4H
W5CHHSX	included	N/A	N/A	W4BAC	included	VR5C	W4BAR	FSM	EA3	PS4H